

Tip
#1

The Fill Line goes on the same side as the DANGER diamond



Since the sharps container box is rectangular it is easy to make the mistake of putting the lid on backwards. When installed correctly,

the Fill Line should be on the same side as the DANGER diamond.

This ensures that the container handle is not an obstruction to the user when inserting medical sharps.



Tip #2

Make sure to secure the lid to the base

This is the
container lid

This is the
container
base



Remember to secure the lid to base of the container, by **pushing down on each corner of the lid until you hear a click.**

Need more help ? Refer to our [linked video](#) to view how to assemble the New Sharps Container.



Tip
#3

How to temporarily close your container



This is the
rectangular
tab



Lightly press the container lid shut until you hear a click. **The rectangular tab will remain in an upright position for you to easily reopen your container.**

Need more help ? Refer to our [linked video](#) to view how to assemble the New Sharps Container.



Tip #4

How to permanently close your container

These are the 2 divots



Press down firmly on the 2 textured divots located on the center of the lid until you hear a click. **When closed permanently, the rectangular tab should no longer be accessible to open the container.**

Need more help ? Refer to our [linked video](#) to view how to assemble the New Sharps Container.

